



Connect Your 3Shape Trios® to DIL Dental Lab's Workflow

Save time and money when you send scans directly from your 3Shape Trios® digital impression scanner.

- In your 3Shape Communicate Account, add cadcam@dildentallab.com to your list of labs.
- Call us at 501-352-6580 to accept your connection. DIL Dental Lab must authorize your connection to our workflow before you can transfer files.
- Once we confirm your request, visit Laboratory Selection and select DIL Dental Lab as your laboratory.

If you have any questions, contact DIL Dental Lab at 501-352-6580 and ask for digital impression technical support.

Looking for more information on our products?

For material descriptions, preparation and cementation guides and billing codes, visit dildentallab.com

Save \$20 on every digital submission by going modeless.

\$79.00

BruxZir®

Full-Contour Zirconia Single Unit



Special introductory price 10/1/16 through 12/31/16

\$99 regular price

DIL Dental Lab • 520 West Pershing, Suite F • North Little Rock, AR 72114 • 501-352-6580 • dildentallab.com

